## ENGLISH - INSTRUCTIONS FOR USE PLEASE READ CAREFULLY MAKING SURE YOU FULLY UNDERSTAND THEM.

Please ensure that you read all the instructions in this leaflet before attempting to wear your new eye accessories. By following these instructions you will get the most from your new look.

Preparing to wear Hygiene is an essential part of preparation. Always wash and dry your hands thoroughly, removing trace of anything that may cause contamination or damage, before handling your eye accessories. Always soak your eye accessories for at least 2 hours before first use in a suitable multipurpose soft contact lens solution.

### Achieving your look

Ensure your eye accessories have soaked prior to first wear. Try to make a habit of putting on the right side first and always make a distinction between right and left. Place one eye accessory on your index finger and check that it is clean and not inside out, (ensure that it has a natural bowl like shape). Hold your upper lid open with your left hand's index and middle finger, at the same time hold your lower lid open with the middle finger of your right hand. Looking straight at your eye in the mirror, slowly and gently place the accessory on your eye. Carefully release your eyelids and blink. You can gently massage your eyelid to centre the accessory. When you are sure that it is correctly inserted, repeat the procedure for the other eye.

### Removing

Wash and dry your hands thoroughly. Look upwards and hold down the lower lid with the middle finger of your right hand. Placing your index finger and thumb of the same hand on opposite sides of the accessory, slowly slide it down to the white of your eye. Gently pinch the accessory from your eye, taking care not to pinch too hard. Repeat for the other eye.

### Wearing instructions

Eye accessories are intended for occasional use only. Check with your eyecare practitioner before wearing for prolonged periods on an everyday basis. Your eyecare practitioner will recommend a suitable wearing schedule. It is suggested that you only wear your eye accessories for 4 hours on the first day, increasing your wearing time by 2 hours daily until you reach optimum wearing time. If you have not worn your eye accessories for several days, you should follow wearing schedule recommendations as a new wearer. Always rinse your accessories, before wearing and after removing, with a suitable multipurpose soft contact lens solution. After each wear, always store your eye accessories in fresh multipurpose soft contact lens solution, making sure not to damage, or catch your eye accessories when storing in any lens case used. If you wear your eye accessories infrequently, always remember to replace the solution used with fresh solution regularly, and ensure that the case used for storage of your accessories is cleaned daily.

These eye accessories have a maximum number of days use (see outer label) and are intended to be replaced at this time after opening, irrespective of the number of times worn.

### Safety

DO NOT wear if your eye accessories appear to be damaged. DO NOT share your eye accessories with anyone else. DO NOT let hairspray or any other cosmetic aerosols come into contact with your eye accessories. DO NOT participate in any water sports whilst wearing. DO NOT wear your eye accessories when using eye medication. DO NOT use tap water, saliva or any household cleaning agents to clean or wet your eye accessories or storage case. ALWAYS remove your eye accessories before you go to bed. ALWAYS put your eye accessories in before applying any make-up. ALWAYS take out eye accessories before removing any make-up. ALWAYS replace your eye accessories if your eyes get red or irritated or if you have any discomfort or pain, and seek help from a qualified eyecare practitioner.

### Important information

Please note that these Eye Accessories are not intended to correct any defect of sight and are intended only for cosmetic colour change.

If you suffer from any of the following, then you MUST NOT wear these eye accessories : Corneal dystrophies Corneal ulcers Inflammation of the cornea, e.g. keratitis Inflammation of the eye, e.g. itis, uveitis Allergic eye conditions Any pre-existing eye condition if you are in any doubt, please seek further advice from medical or eyecare practitioner. Eye Accessories are suitable for normal, healthy eyes. If you are in any doubt as to whether these eye accessories are suitable for you, then please consult a medical or eyecare practitioner.

### Fitting and Aftercare

It is recommended that you consult your eyecare practitioner to ensure suitability and for practical advice on wearing and removing.(consultation fees may be charged).

Always have regular check-ups and follow the advice of your medical or eye care practitioner. If you feel any unexplained discomfort, watering, vision change or redness,remove immediately and seek help from a qualified eyecare practitioner.

If you have any questions or concerns regarding the information in this leaflet, Please consult your medical or eyecare practitioner.

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